

New vmed Vet-Tab

Miniature Wireless Patient Monitor

Bluetooth to Android Device

Our least expensive wireless monitor has built-in display and connects to your Android device for remote viewing.



\$1,995 with free tablet

Validated for small animals

- Noninvasive blood pressure
- ECG three lead
- SpO2 with lingual clips
- Temperature Celsius and Fahrenheit
- Respiration from ECG

Features

- Small size, weighs 9 ounces
- Built-in color display
- Monitor and control from Android app
- Easy one button operation
- Time between charges 12 hours
- Mounting hooks for surgical table and kennel
- Settings menu
- Optional esophageal probes and other sensors

Warranty

Two years parts and labor
Six months for accessories

Battery and AC adapter

Battery type: 3.7V,1800mAh, lithium-ion.
Time between charges 9-12 hours
AC adapter: Input: AC 100-240v
~200mA, Output: DC 5V 1A

Parameters Specifications

ECG Three Lead

Band width: diagnostic mode 0.05Hz-100Hz (3dB), surgery mode: 0.1-25Hz
Heart rate range 15 bpm to 300 bpm
Accuracy: ± 1 bpm

SpO2

Detection range: 35%-100%
Accuracy: $\pm 1\%$ 90%-100%, $\pm 2\%$ 70%-89%

NIBP


Cuff inflation time: <15 seconds
Initial inflation pressure: 180 mmHg
Detection range
Systolic 30mmHg-255mmHg
Diastolic 15mmHg -220mmHg
MAP 20mmHg-235mmHg

Respiration

Chest impedance method
Range 0 bpm-100 bpm
Accuracy ± 1 bpm

Temp

Type: Thermistor
Range: 20C-50C
Accuracy: ± 0.2 C



Integral display with one button controls at bottom and indications at top.

vmed technology

vmedtechnology.com, johnb@vmedtech.com
800•926•9622