



MSS McALLISTER
Software Systems®

COMPATIBILITY
PARTNER



Vmed Technology

Vmed Technology Integrations within AVImark:

The Vmed integration allows for specific client/patient data to be pulled from the AVImark system and automatically populated in Vmed's "patient episode file" that is created by the PC Display software. This software is used in conjunction with Vmed's PC-VetGard and PC-Vet wireless monitors for the purpose of storing patient monitor parameters and ECG waveform data.

Vmed Technology

- PC Display
- PC-VetGard & PC-Vet

AVImark version

version 2009.0 & higher

-
- AVImark technical support is available at 877-838-9273.
 - Users may log into on-line chat support at avimark.net
 - AVImark is the world leader in veterinary practice management software with over 9,000 clinics and 20 years of experience.

